

SIZING SYSTEM

PERFECT-FIT*

VenoTrain micro

VenoTrain soft / soft S

		XS for VenoTrain soft S		S		M		L		XL	
				NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
		Übermööß									
		cT		up to 85		up to 90		up to 100		up to 120	
		cH		up to 95		up to 105		up to 120		up to 140	
		Silicone top band, wide**		51-57	59-66	56-63	62-69	60-67	66-74	63-71	71-80
		cG Silicone top band, regular		44-50	52-59	49-56	55-62	53-60	58-66	55-63	62-71
		Top band, AT, silicone top band		44-50	52-66	49-63	55-69	53-67	58-74	55-71	62-80
		cF		40-52	43-55	43-55	46-60	45-59	50-64	49-63	52-68
		cE		33-40	36-44	35-43	39-49	38-48	42-53	42-53	47-58
		cD		29-35	32-38	32-38	35-42	34-42	37-45	37-45	41-49
		cC		31-37	34-40	34-40	37-44	37-45	39-47	40-48	43-51
		cB1		25-29	27-31	27-32	29-34	29-35	31-37	32-38	34-40
		cB		20-22		22-25		25-28		28-31	
		cY		29-33		31-35		34-40		37-44	
		cA		19-24		22-27		25-31		28-35	

Length lG

short 61-70

long 70-80

Type AG, AT

Length lD

short 35-40

long 40-46

Type AD